

George C. Marshall Space Flight Center
Marshall Space Flight Center, Alabama 35812

ED27-EMA-FOP-010 Revision B August 28, 2001

FACILITY OPERATING PROCEDURE

ED27 / MODAL AND CONTROL DYNAMICS TEAM

VERIFICATION OF DATA ACQUISITION SYSTEM USED IN MODAL TEST LAB

CHECK THE MASTER LIST—
VERIFY THAT THIS IS THE CORRECT VERSION BEFORE USE

ED27 / Modal and Control Dynamics Team			
Verification of Data	ED27-EMA-FOP-010	Revision: B	
Acquisition System Used in			
Modal Test Lab	Date: 8/28/01	Page 2 of 5	

DOCUMENT HISTORY LOG

Status (Baseline/			
Revision/	Document	Effective	
Canceled)	Revision	Date	Description
Baseline		7/28/99	Document converted from ED73-EMA-FOP-010 Revision B due to reorganization. Previous history retained in system as part of canceled Group files. Revised entire document to reflect the new organizational structure and document numbering system.
Revision	A	8/8/00	Update to include new DIFA Scadas III
Revision	В	8/28/01	Update to include HP VXI

ED27 / Modal and Control Dynamics Team			
Verification of Data	ED27-EMA-FOP-010	Revision: B	
Acquisition System Used in			
Modal Test Lab	Date: 8/28/01	Page 3 of 5	

1. INTRODUCTION

- 1.1 <u>Scope</u> This procedure defines the system used to verify the data acquisition systems used in the performance of experimental modal analysis and testing.
- 1.2 <u>Purpose</u> This procedure defines the system to fulfill the requirements of $\overline{\text{ED27-OWI-M\&V-002}}$ "Quality Records Control".
- 1.3 <u>Applicability</u> This procedure applies the data acquisition systems used in the performance of experimental modal analysis and testing within the Modal and Control Dynamics Team/ED27.
- 1.4 <u>Frequency</u> This procedure shall be performed on all Data Acquisition Systems (Section 3.1) at one (1) year intervals.

2. REFERENCE DOCUMENTS

_
the HP9000 Computers with
dware for Modal Surveys.
the HP9000 Computers with
ement Hardware for Modal
the NT Workstation Computer
surement Hardware for Modal

3. DEFINITIONS

- 3.1 <u>Data Acquisition System</u> The test equipment known as the "Data Acquisition System" is composed of the following components.
 - 3.1.1 <u>Signal Conditioners</u> Consists of one or more instruments which provide a power source to the transducer and provide amplification of the output signal.
 - 3.1.2 <u>Measurement Hardware</u> Consists of one or more modules which provide analog-to-digital conversion of the transducer output signal and sends the digital data to the host computer.
 - 3.1.3 <u>Host Computer</u> Provides a control interface to the measurement hardware and allows for manipulation and data storage of the acquired data.
- $3.2 \ \underline{\text{Auxiliary Equipment}}$ Describes any other equipment used in the verification process that is not part of the Data Acquisition System and includes the following components.
 - 3.2.1 Hand-Held Calibration Exciter w/Reference Accelerometer
 Provides a constant amplitude and frequency signal as input to the
 Data Acquisition System. The Hand-Held Calibration Exciter is
 considered to be Category I test equipment (Mandatory Calibration
 and Recall) and its calibration is documented as such. The
 Reference Accelerometer is a Category II item.

ED27 / Modal and Control Dynamics Team			
Verification of Data Acquisition System Used in	ED27-EMA-FOP-010	Revision: B	
Modal Test Lab	Date: 8/28/01	Page 4 of 5	

- 3.2.2 <u>Cabling and Patch Panel</u> Provides an interface between the Reference Accelerometer and the Signal Conditioner. Both of these items are Category III items (Not Calibrated).
- 3.2.3 <u>Digital Multimeter</u> Measures both the amplitude and frequency of the reference signal. The multimeter is considered to be Category I test equipment (Mandatory Calibration and Recall) and its calibration is documented as such.

4. INSTRUCTIONS

- 4.1 <u>Preparation</u> Set up and cable the Data Acquisition System as described in the applicable reference document.
- 4.2 <u>Verification Procedure</u> This procedure describes the method used to verify the data acquisition system.
 - 4.2.1 Cable the Reference Accelerometer or reference signal to the first channel of the Data Acquisition System.
 - 4.2.2 Initiate the data acquisition process from the host computer.
 - 4.2.2.1 From the FMON/MIMO select Channels..., Identification... and set the channel dimension to "Voltage".
 - 4.2.2.2 Then select Channels..., Calibration..., Absolute.
 - 4.2.2.3 Enter the expected output frequency from the hand-held calibrator or reference signal in the "Reference" field.
 - 4.2.2.4 Type the expected voltage output by the reference accelerometer or reference signal in the "Amplitude" field.
 - 4.2.2.5 Set the Reference Value Type to "RMS".
 - 4.2.2.6 Choose the Autorange and Calibrate option.
 - 4.2.2.7 Select desired channel (by "rubber banding") and hit the "Start" button.
 - 4.2.2.8 After the calibration cycle is complete, read the value in the "Calibration" field and transfer to quality record.
 - 4.2.3 Repeat for all channels of the Data Acquisition System.
 - 4.2.4 If the percent variance for frequency exceeds %2.0 or the percent variance for amplitude exceeds 4.0% for a given channel, the input module containing that channel will be taken out of service and repaired.

5. QUALITY RECORDS

5.1 <u>Calibration Log Book</u> The following information will be recorded in the calibration log book as an electronic database per ED27-OWI-M&V-002: Model number and NASA property number of the digital multimeter, date this procedure was performed, the measured amplitude and frequency of each channel as recorded in section 4.2.2 above, and the percent variance for both the recorded frequency and amplitude.

ED27 / Modal and Control Dynamics Team			
Verification of Data	ED27-EMA-FOP-010	Revision: B	
Acquisition System Used in			
Modal Test Lab	Date: 8/28/01	Page 5 of 5	

Table A: Sample Calibration Record

Ch #	Amp (g's)	Amp % var.	Freq (Hz)	Freq % var.	Date
1	0.997	0.2%	79.625	0.15%	3/6/98
2	0.996	0.3%	79.625	0.15%	3/6/98
•••		•••	•••	•••	
16	0.993	0.6%	79.625	0.15%	3/6/98

Keithley 2000 M628880 due: 8/27/00